

Redefining procedural efficiency


without compromising tissue quality - from positioning to diagnosis.



Mammotome Revolve™ dual vacuum-assisted breast biopsy system for stereotactic-guided procedures provides efficiency and clinical confidence throughout the procedure and onto pathology. Its advanced tissue management system automatically organizes specimens into touch-free collection chambers, which are designed for ease of specimen x-ray and transport to pathology while preserving tissue integrity. Patented dual vacuum technology acquires the desired volume of tissue with fewer passes, while maximizing the amount of intact tissue for clinical confidence.









Provides touch-free tissue collection into numbered, pathology-ready chambers. Cartridge trays can be flattened and imaged to identify and easily locate calcifications. Organized presentation and visualization of the tissue specimens allow real-time sampling adjustment. Confirm sample acquisition at any needle position.





Dual Vacuum Technology

Continually holds and secures tissue within

aperture resulting in larger, contiguous samples

delivering a more complete picture to pathology.

Optimal Firing Speed

Firing speed is designed to minimize lesion

displacement and increase targeting accuracy.




Redefining procedural efficiency

without compromising tissue quality - from positioning to diagnosis.



Mammotome Revolve™ dual vacuum-assisted breast biopsy system for ultrasound-guided procedures provides efficiency and clinical confidence throughout the procedure and onto pathology. Its robust, yet lightweight and ergonomic handpiece fits comfortably in your hand, and its advanced needle technology improves ease of insertion. Patented dual vacuum technology acquires the desired volume of tissue with fewer passes, while maximizing the amount of intact tissue for clinical confidence.





Specimen Management System

Provides touch-free tissue collection into numbered, pathology-ready chambers. Organized presentation and visualization of the tissue specimens allow real time sampling adjustment. Confirm sample acquisition at any needle position.







Designed to empower clinical confidence

Order Placement

Capital Equipment

| Code | Description | Quantity |
|---------|--|----------|
| MSCM1 | Mammotome Revolve™ System Control Module | 1 |
| MCART1 | Mammotome Revolve™ System Cart | 1 |
| MHKEYP1 | Mammotome Revolve™ Remote Keypad | 1 |
| MF00T1 | Mammotome Revolve™ Remote Footswitch | 1 |
| MSTH1 | Mammotome Revolve™ Stereotactic Holster | 1 |
| MHUSH1 | Mammotome Revolve™ Ultrasound Holster | 1 |

Mammotome® Markers

| Code | Description | Shape | Quantity |
|---------------|--|------------------|----------|
| 4010-05-08-T1 | 8G HydroMARK™ Breast Biopsy Site Marker | | 10 |
| 4010-05-08-T3 | 8G HydroMARK™ Breast Biopsy Site Marker | INN | 10 |
| 4010-05-08-T4 | 8G HydroMARK™ Breast Biopsy Site Marker | | 10 |
| 4010-05-10-T1 | 10G HydroMARK™ Breast Biopsy Site Marker | | 10 |
| 4010-05-10-T3 | 10G HydroMARK™ Breast Biopsy Site Marker | INN | 10 |
| 4010-05-10-T4 | 10G HydroMARK™ Breast Biopsy Site Marker | | 10 |
| MMK0801 | 8G MammoMARK®/CorMARK® Biopsy Site Identifier | \triangleright | 10 |
| MMK0802 | 8G MammoMARK®/CorMARK® Biopsy Site Identifier | | 10 |
| MMK0803 | 8G MammoMARK®/CorMARK® Biopsy Site Identifier | 0,00 | 10 |
| MMK1001 | 10G MammoMARK® / CorMARK® Biopsy Site Identifier | | 10 |
| MMK1002 | 10G MammoMARK® / CorMARK® Biopsy Site Identifier | | 10 |
| MMK1003 | 10G MammoMARK® / CorMARK® Biopsy Site Identifier | 0,00 | 10 |
| STAR0831 | 8G MammoSTAR® Biopsy Site Marker | \subset | 10 |
| STAR0832 | 8G MammoSTAR® Biopsy Site Marker | \sim | 10 |
| STAR0833 | 8G (P) MammoSTAR® Biopsy Site Marker | \subset | 10 |
| STAR1031 | 10G MammoSTAR® Biopsy Site Marker | \subset | 10 |
| STAR1032 | 10G MammoSTAR® Biopsy Site Marker | \iff | 10 |
| STAR1033 | 10G (P) MammoSTAR® Biopsy Site Marker | \approx | 10 |

Stereotactic Probes

| | Code | Description | Quantity |
|--|---------|--|----------|
| | MST0809 | 8G x 9cm Mammotome Revolve™ Stereotactic Probe | 5 |
| | MST0812 | 8G x 12cm Mammotome Revolve™ Stereotactic Probe | 5 |
| | MST0815 | 8G x 15cm Mammotome Revolve™ Stereotactic Probe | 5 |
| | MST1009 | 10G x 9cm Mammotome Revolve™ Stereotactic Probe | 5 |
| | MST1012 | 10G x 12cm Mammotome Revolve™ Stereotactic Probe | 5 |
| | MST1015 | 10G x 15cm Mammotome Revolve™ Stereotactic Probe | 5 |
| | | | |

^{*}Probe includes Tube Set and Specimen Management System

Ultrasound Probes

| MHUS08 | 8G Mammotome Revolve™ Ultrasound Probe | 5 |
|--------|---|---|
| MHUS10 | 10G Mammotome Revolve™ Ultrasound Probe | 5 |

Send your order via fax to your Customer Service Representative at

Disposable Accessories

| Code | Description | Quantity |
|------------|---|----------|
| MG08A | 8G Mammotome Revolve™ Probe Guide, Fork Shape | 25 |
| MG10A | 10G Mammotome Revolve™ Probe Guide, Fork Shape | 25 |
| MG08B | 8G Mammotome Revolve™ Probe Guide, Button Shape (for GE Upright Systems) | 25 |
| MG10B | 10G Mammotome Revolve™ Probe Guide, Button Shape (for GE Upright Systems) | 25 |
| MCANISTER1 | Replacement Vacuum Canister (800cc) for Mammotome Revolve™ System | 10 |
| MSMB1208 | 8G (12 chamber) Replacement Specimen Management System | 10 |
| MSMB1210 | 10G (12 chamber) Replacement Specimen Management System | 10 |







^{*}Probe includes Tube Set and Specimen Management System

Differentiation at Each Step of the Stereo Process



Tissue AcquisitionEliminate tissue handling



Specimen Imaging/ConfirmationConfirm at the point of care



Biopsy Site MarkingEnable post-biopsy site identification

- Reduce procedure time
- Enhance patient experience
- Improve resource utilization

Mammotome

mammotome.com

For complete product details, see Instructions for Use.



Product may not be approved or available in your region. Please check with your local Mammotome representative.