

Advancing Cancer Diagnostics
Improving Lives

Leica
BIOSYSTEMS

The Compact Solution to Quality Tissue Processing



**HistoCore PEARL
Tissue Processor**

HistoCore PEARL Tissue Processor




Ultimate Protection for User Safety

- Fully enclosed system
- Triple Fume Protection: extraction at the retort, an external exhaust adapter and an activated carbon filter provide added safety for lab staff



Ultra-Compact Solution

- 200 cassette capacity
- Small footprint to maximize laboratory space
- 3.5L reagent bottles and three paraffin baths help maximize run efficiency



Simple Operation

- Easy-to-remember operating procedures on xylene resistant, multi-lingual touch-screen interface
- Three pre-validated protocols and Reagent Management System (RMS) provides added flexibility
- USB connection to download log files help maintain accurate records



Easy-to-Clean Paraffin and Reagent Exchanges

- Drawer style design for paraffin baths allows for easy filling and draining
- Trays under the paraffin baths and reagent bottles effectively collect spilled liquids/chemicals





Contact your Leica Biosystems representative for a full range of validated consumables for the HistoCore PEARL Tissue Processor.

SYSTEM SPECIFICATIONS

GENERAL

Nominal supply voltages	100 to 120 V AC or 220 to 240 V AC
Nominal frequency	50 / 60 Hz
Nominal power	1650 VA
Operating temperature range	+10 °C to +35 °C

RETORT

Capacity	Max. 200 cassettes
Temperature (paraffin)	50 °C to 70 °C
Temperature (processing reagents)	Ambient temperature or 35 °C to 60 °C
Temperature (cleaning reagents)	Alcohol 62 °C, Xylene 67 °C

PARAFFIN BATHS

Number of paraffin baths	3
--------------------------	---

Capacity	Max. 3.5 L per paraffin bath
Melting time	50 °C to 64 °C, 12 hours
Temperature setting range	50 °C to 70 °C

BOTTLES

Number of reagent bottles	11 (in reagent cabinet)
Condensate bottle	1
Cleaning solution bottles	2
Bottle volume	Max. 3.5 L

DIMENSIONS AND WEIGHTS

Dimensions, unpacked (W x D x H)	600 x 710 x 1320 mm (23.62 x 27.95 x 51.96 in)
Weight, unpacked	Max. 140 kg (308.65 lbs)
Weight, packed	Max. 175 kg (385.81 lbs)

Leica Biosystems—an international company with a strong network of worldwide customer service.

North America Sales and Customer Support

North America	800 248 0123
---------------	--------------

Asia/Pacific Sales and Customer Support

Australia	1800 625 286
China	+85 2 2564 6699
Japan	+81 03 6758 5690
South Korea	+82 2 3416 4500
New Zealand	0800 400 589
Singapore	+65 6550 5999

Latin America Sales and Customer Support

Brazil	+55 11 2764 2411
--------	------------------

Europe Sales and Customer Support

For detailed contact information about European sales offices or distributors, please visit our website: LeicaBiosystems.com

Contact your Leica Biosystems representative today to learn more about the HistoCore PEARL Tissue Processor.

LEICA BIOSYSTEMS

Leica Biosystems is a global leader in workflow solutions and automation, integrating each step in the workflow from biopsy to diagnosis.

**Advancing
Cancer Diagnostics,
Improving Lives**

